

Materials In Medicine

Advances in Gowns & Gloves



FELIX

Eskisehir Industrial Zone
26th Street No:9
26101 , Eskisehir - Turkey

Phone : +90 222 236 23 63
Fax : +90 222 236 23 53


www.felixplastic.com

felix@felixplastic.com

Technologies for PPE – State of the Art

Material Properties

- Reusable/Disposable
- Treatments/Barriers
- Absorption/Strength




What is necessary?

- A simple and cost effective portfolio of apparel and protective products that is proven to provide the *highest level of protection and comfort* whilst offering *exceptional ease of use*
- Key to success for protection is *practice*, namely:
 - Standard Operating Procedures (SOPs) in place
 - Routine use of PPE
 - Right PPE for the job
 - Wear length of PPE
 - Disposal of PPE

Technologies for PPE – State of the Art

Key PPE necessitates HIGH PERFORMANCE MATERIALS as seen in:

- 1 **Bodygard® Surgical Drapes**
 - Fully impervious surgical drape
 - Absorbent top layer for fluid management
 - AAMI Level 4*
- 2 **Spunsper -SFS™ Surgical Gown**
 - Breathable, fully reinforced surgical gown
 - AAMI Level 4* F1670/F1671
 - 180 cm hydrostatic pressure
- 3 **SFS Bodygard™ Surgical Gown**
 - Breathable fully reinforced comfortable surgical gown
 - AAMI Level 3* 180cm hydrostatic pressure
- 4 **SMS Surgical Gown**
 - AAMI Level 1 and 2, (task-dependent)
 - AAMI Level 3 is available with upto 55cm hydrostatic pressure




AAMI Standard PB70:2003 (Association for the Advancement of Medical Instrumentation)

** Liquid Barrier Performance Classification of Protective Apparel & Drapes Intended for Use in Health Care Facilities*


FELIX

Technologies for PPE – State of the Art


- AAMI defines critical zones as those areas where direct contact is likely to occur with:
 - Blood
 - Body Fluids
 - Other potentially infectious material (OPIM) including influenza virus - e.g., H1N1, HCV
- Critical areas include material as well as areas of construction - e.g., seams



Drape Reinforcement




Isolation Gown / Fully
reinforced surgical
gown(entire gown)



Critical Zones

Surgical Gown

Technologies for PPE – State of the Art




- **Three-layer micro-fiber composite:**
 - Absorbent top layer
 - Impervious middle membrane
 - Patient comfort layer
- **Key attributes:**
 - Impervious (AAMI Level 4)
 - Flame-resistant
 - Absorbent
 - Abrasion resistant and low lint
 - Puncture resistant



1-Bodygard® Surgical Drapes

FELIX

Technologies for PPE – State of the Art



Water Vapor Transmission Rate (WVTR)



Superior protection against bacteria & viruses (H1N1, HCV)

2-Product Performance - Spunsper-SFS Bodygard™

FELIX

Technologies for PPE – State of the Art

3-SFS *Bodygard*[™] Surgical Gowns


- Soft, silky material makes it a comfortable gown
- AATC 127 upto 180cm hydrostatic pressure
- AAMI Level 3 protection
- ***60gsm Lighter material helps reduce surgical glove roll down***



FELIX

Technologies for PPE – State of the Art


Isolation Gowns



- Frequently used form of PPE
- Requires FDA* approval if it provides enhanced barrier protection – e.g., antimicrobial, antiviral
- Disposable – single-use
- Intended for short-term use or procedures when there is a risk of exposure to blood, body fluids, OPIM (e.g., H1N1), chemicals, or chemotherapy drugs
- Four common types:
 - PE Film
 - Spunbond
 - SMS
 - SFS (Spunbond+film+Spunbond)
- 2004 – 2009: Year-on-Year Growth is 5%

* Draft Guidance for Industry and FDA Staff: Premarket Notification (510K) submissions for medical devices that include antimicrobial agents – July 19, 2007

Technologies for PPE – Current & Future

Gown Type	Usage	PROS	CONS	FUTURE
Spunbond - most appropriate for low-risk procedures that require limited barrier protection.	<ul style="list-style-type: none"> •Visitors •Zero fluid contact •Used in various wards •Patient care > 5+ min 	<ul style="list-style-type: none"> •High comfort and durability •Inexpensive 	<ul style="list-style-type: none"> •Single-use •Limited protection 	<ul style="list-style-type: none"> •Very basic educate on protection to move users to SMS based gowns
SMS appropriate for mid-risk procedures that require some additional barrier protection.	<ul style="list-style-type: none"> •Changing bed linens •Nursing/Administering an injection •Used in ICU's /NICU's/PACU's •AAMI Level 1,2, or 3 •Patient care > 15+ min 	<ul style="list-style-type: none"> •Improved comfort and durability •Medium Cost 	<ul style="list-style-type: none"> •Single-use •Mid-Level protection 	<ul style="list-style-type: none"> •Include additional barrier protection to impact frequency/transmission of HAIs, OPIM (H1N1) & MRDOs •'Green' for disposal
PE Film - used for procedures that require high levels of barrier protection.	<ul style="list-style-type: none"> •Janitors/Auxiliaries •One time use •Sterile Processing Dept •Patient care < 15 min 	<ul style="list-style-type: none"> •Inexpensive •High Barrier 	<ul style="list-style-type: none"> •Non-breathable •Uncomfortable •Flimsy 	<ul style="list-style-type: none"> •Very basic – move to SMS +PE gown
SFS - best for procedures requiring a high level of barrier protection	<ul style="list-style-type: none"> •Preparation & administrations of chemotherapeutic agents •AAMI level 4 •Patient care > 15+ min 	<ul style="list-style-type: none"> •Excellent barrier protection •Good chemical resistance 	<ul style="list-style-type: none"> •Costly 	<ul style="list-style-type: none"> •Barrier improvements for aggressive agents and additional chemodrugs •'Green' for disposal


Technologies for PPE – Evidence Basis for PPE

HAI's - Isolation Gowns & Gloves

- 1.7M acquire & some 100,000 die from nosocomial infections each year*
- Bacteria responsible for spread of such infections: Staphylococci & E-coli.
- Pathogenic agents transmitted via the skin (hands), through droplets (from mouth and throat) and, last but not least, **by textiles**.
- Critical for spread are nurses' and orderlies' tunics and physicians' coats (especially the coat pockets). (*AJIC*, Aug 07, vol 35)
- End-Users (Patient Providers) - points to value of antimicrobial agents as a defense against HAI's & antiviral coatings against H1N1

VRE - Isolation Gowns & Gloves

- Recent study showed routine use of isolation gowns & gloves prevented more than 58 VRE cases during an 18 month period**
- Cost of treating a patient with infection is some \$60K dollars
- Cost of providing gowns/gloves to prevent is approx. \$2K per patient



* **Estimating HAIs and Deaths in US Hospitals, 2002** – Kleven, R.M. et al – Public health Reports, 2007, Vol. 122

** **A cost benefit analysis of Gown use in controlling VRE transmission; is it worth the price?** - Puzniak, L.A., Gillespie, K.N., Leet, T., Kollef, M. and Mundy, L.M. – Infection Control & Hospital Epidemiology, 2004, Vol 25, no.5, p 418-424

Design & Engineering Breakthroughs on Horizon

Materials

- Non-wovens for disposable products
 - with new additives – lighter, more absorbent
 - ‘green’ sustainable – PLA, corn, non-petroleum based resin, algae-based resin
 - with barriers – hydrophobic, antimicrobial, antiviral, radiation, self-sealing, repellency
- Alternative fibers
 - Bamboo
 - Carbon nanotubes (SWCNT)
 - Other cellulosic

Products

- System solutions – e.g., combination device
- Advocate for single-use, disposable
- Trends match automation
- Improved comfort – softer, more breathable
- Easy to don/doff - minimize contamination

Other Technologies

- Indicators
- Sensors
- Smart, Smarter, Smartest



Design & Engineering Breakthroughs on Horizon

Time Zone	Performance Materials	Products Protocols	Systems Practice
1-3 years 'SMART'	<ul style="list-style-type: none">•Lower basis weight•Move to bi-laminates•Odor elimination•Maintain/lower cost	<ul style="list-style-type: none">•Improved fit, comfort•Enhanced barriers - antimicrobial, antiviral•Minimize cross-contamination•Maintain/lower cost	<ul style="list-style-type: none">•Connectivity - gowns to gloves•Support compliance for effective protection
3-5 years 'SMARTER'	<ul style="list-style-type: none">•Move to 'green' materials•Improved additive technology•Monitoring capabilities - temperature, blood pressure and other vital signs•Maintain/lower cost	<ul style="list-style-type: none">•Include more sensor technologies, feedback and alerts•Tear-proof/puncture resistant gown and gloves	<ul style="list-style-type: none">•Easy to don/doff•Dual purpose products - gown & wipe•Combo products - Gown + Glove•Support compliance via RFID tracking systems
5-10 years 'SMARTEST'	<ul style="list-style-type: none">•Fibers and shapes for performance vs. cost•Functionalization of surfaces – electrical and chemical activity and responsiveness	<ul style="list-style-type: none">•Self-cleaning•Maintain/lower cost	<ul style="list-style-type: none">•New technologies for improved cleaning in hospitals•Reduced bio-burden on PPE for infection control

Ultimate Goal: Make PPE Obsolete



Thank you!